


# Reusable SpO<sub>2</sub> Sensor

## Instructions for Use

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## English SpO<sub>2</sub> SENSOR INSTRUCTIONS FOR USE

#### Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this product and this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observance of this manual is a prerequisite for proper performance and correct operation, and ensures patient and operator safety.

- Mindray is responsible for safety, reliability and performance of this product only in the condition that:
  - This product is operated under strict observance of this manual.
  - This product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.
- In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number must appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief description of the reason for return. The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

**mindray**  **MINDRAY** and **A** are the registered trademarks or trademarks owned by Mindray in China. The issue date of this manual is 2021-06. © 2007-2021 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

**Product Introduction**
**Intended Use**
The reusable 518B, 518BL and 518BLH SpO<sub>2</sub> sensors are intended for noninvasively spotchecking or continuously monitoring the arterial oxygen saturation and pulse rate of neonatal patients (weight lower than 5 kg).
The 518B, 518BL and 518BLH SpO<sub>2</sub> sensors must be used in connection with the monitoring device (including accessories) specified by our company, including the MEC, IMEC, PM, VS, and T series patient monitors of Mindray that supports the SpO<sub>2</sub> measurement function. The user shall be responsible for verifying the compatibility of this sensor, the extension cable, and the equipment before use. This sensor must be used by trained clinical professionals. Otherwise, product performance may be compromised or personal injury may occur.
This sensor material that contacts the patients has undertaken the bio-compatibility test and is verified to be in compliance with ISO 10993-1.
Select the corresponding SpO<sub>2</sub> extension cable based on the type of SpO<sub>2</sub> connector.

#### Contraindications


None.

#### Warnings

- Federal law (USA)** restricts this device to sale by or on the order of a physician.
- Check the finger sensor and the cable before using. Do not use them for SpO<sub>2</sub> monitoring in case of damage, block, deterioration or visible contamination.
- For disposal of accessories, follow local regulations or your hospital's policy regarding disposal of such accessories. Do not dispose randomly.**
- Where possible, apply the sensor to the site free of arterial catheters, blood pressure cuffs, or intravascular infusion lines.
- Interference of high frequency electricity may affect the accuracy of the measurement.
- Try to keep the patient still (especially the arm) and avoid the measured site suffering excessive motion.**
- Cover the sensor with opaque material under the condition of strong light. Failure to do so will result in inaccurate measurement.**
- Always choose a site that is well perfused and will completely cover the sensor's detector window. Site should be cleaned of debris and dry prior to sensor placement.**
- The SpO<sub>2</sub> sensor must be used in connection with the specified equipment only. The user shall be responsible for ensuring the compatibility between this sensor, the extension cable and other equipment by either referring to the instructions for using this sensor or contacting us before use. Otherwise, product performance may be compromised or personal injury may occur.**
- Carefully route cables to reduce the possibility of patient entanglement or strangulation.**
- Do not use SpO<sub>2</sub> sensor during MRI scanning. Conducted current may cause burns.**
- The sensor application site must be inspected every two hours.**
- With very low perfusion at the monitored site, the reading of oxygen saturation may be inaccurate.**
- If using SpO<sub>2</sub> sensor during full body irradiation, keep the sensor out of the irradiation field. If the sensor is exposed to the irradiation, the reading might be inaccurate or the product might reach zero for the duration of the active irradiation period.**
- Always select the site according to different patients and application sites. Tight sensor may cause venous pulsation, obstructed blood circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurement, while loose sensor may lead to erroneous optical alignment or falling off. If the sensor is too tight because the application site is too large or becomes too large due to edema, excessive pressure may result in congestion distal from the application site, leading to interstitial edema and tissue ischemia.**
- Avoid bending the joint of the sensor cable for a long time. Otherwise, it will reduce the life of sensor cable.**
- Do not place the product in an environment that exceeds the stated range.**
- Do not bend the cable excessively for a prolonged time.**
- When performing continuous long-term monitoring on a patient, regularly inspect and change the sensor application site according to the skin quality. Inspecting at least every two hours is recommended.**

#### Applying SpO<sub>2</sub> Sensor

- Be sure to understand all warnings listed in the previous section before applying any sensor to a patient.
- Make sure the sensor is free of dirt and rust before use. Clean the sensor or replace it if necessary. Check the compatibility between the sensor and the monitor. Incompatibility could affect the measurement accuracy, performance, and reliability of the product. It may even bring harm to the patient. Apply the SpO<sub>2</sub> sensor to the patient, as shown below.
  - For 518B and 518BL SpO<sub>2</sub> sensors, clip the neonatal foot with the sensor. Pull the belt and allow it to go through the corresponding belt fastener and then through the fastener on the other side of the sensor jacket to fix the SpO<sub>2</sub> sensor.
  - For 518BLH SpO<sub>2</sub> sensor, clip the neonatal heel with the sensor. Use the magic tape on the back of the belt to fix the SpO<sub>2</sub> sensor, and then pull the belt and allow it to go through the belt fastener.

Connect the sensor to the cable. Plug the cable into the SpO<sub>2</sub> connector of an instrument in the direction of the blue arrow on the cable.

#### Note:

- Improper sensor placement may affect the measurement accuracy.**
- The sensor should be placed on the same site and the heart rate at the same level.**
- Applying the sensor too tightly may cause inaccurate measurement readings.**
- Ensure that the emitter aligns the photodetector directly and all light generated by the emitter passes through the tissues.**
- Functional tester or patient simulator shall not be used to validate the SpO<sub>2</sub> accuracy of pulse oximeter equipment.**

**Checking the Functionality of the Sensor**
The sensor is abnormal and needs to be changed if the following situations occur when the sensor is connected to the instrument:

- 1) The "Sensor Off" alarm occurs.
- 2) The "Senior Low" OI<sub>2</sub> alarm when the sensor is closed but not applied to the finger.
- 3) The "Sensor OI<sub>2</sub>" alarm when the sensor is applied to the finger.

Start SpO<sub>2</sub> monitoring only when you are sure the sensor functions normally.

- Factors affecting measurement accuracy**
  - The measurements also depend on absorption of special wavelength ray by oxidized hemoglobin and deoxyhemoglobin. Concentration of nonfunctional hemoglobin may affect the accuracy of the measurement.
  - Shock, anemia, hypothermia and the application of vasoconstriction drug may decrease arteria blood flow to an unmeasurable level.
  - Pigment, or deep color (for example: nail polish, artificial nails, dye or pigmented cream) may cause inaccurate measurements.
- Cleaning/Disinfection**
  - **Do not immerse the connector of 518B and 518BLH SpO<sub>2</sub> sensor in any cleaning solutions, disinfectants, or other liquids. The other parts of the sensor can be immersed in the recommended cleaning solutions and disinfectants.** After it becomes moistened (immersed for a maximum of thirty seconds), take it out of liquid immediately and let it air dry.
  - **Exercise caution during cleaning/disinfection to avoid wetting the pins.**
  - **Do not clean or disinfect the metal parts of this product. Otherwise, metal part corrosion or product reliability deterioration could result.**
  - **We recommend that the sensor be disinfected only when necessary as determined by hospital infection control procedures.**
  - **Never use cleaning agents/disinfectants other than the recommended.**





#### Cleaning/Disinfection

- Warning**
  - **Do not immerse the connector of 518B and 518BLH SpO<sub>2</sub> sensor in any cleaning solutions, disinfectants, or other liquids. The other parts of the sensor can be immersed in the recommended cleaning solutions and disinfectants.** After it becomes moistened (immersed for a maximum of thirty seconds), take it out of liquid immediately and let it air dry.
  - **Exercise caution during cleaning/disinfection to avoid wetting the pins.**
  - **Do not clean or disinfect the metal parts of this product. Otherwise, metal part corrosion or product reliability deterioration could result.**
  - **We recommend that the sensor be disinfected only when necessary as determined by hospital infection control procedures.**
  - **Never use cleaning agents/disinfectants other than the recommended.**
- Cleaning**
  1. Clean the sensor with cotton or soft cloth moistened with water. The 518B and 518BLH SpO<sub>2</sub> sensors can also be rinsed and wiped in water.
  2. After cleaning, wipe off the water with a soft cloth.
  3. Allow the sensor to air dry.
- Disinfection**
For recommended disinfectants, see the agents with  in the 512/S12/S18/518BL/518BLH/51 3A/51B column of the *Applicable Disinfecting Agents List* (P/N: 046-011532-00).

  1. Clean the sensor as instructed above.
  2. Immerse the sensor in cotton or soft cloth moistened with one of the recommended disinfectants. Except for the connector, other parts of the 518B SpO<sub>2</sub> sensor and 518BLH SpO<sub>2</sub> sensor can be immersed in recommended disinfectants. For 518BLH SpO<sub>2</sub> sensor, only the belt part with the magic tape can be immersed in recommended disinfectants.
  3. After disinfection, be sure to wipe off the disinfectant left on the sensor with a soft cloth.
  4. Allow the sensor to air dry.

#### SPEC/Symbol

The length of the sensor cable is 1.2 m. The length is approximately 3 m after the cable gets connected with the SpO<sub>2</sub> sensor.
Degree of protection against harmful ingress of water: IPX7 (Protection against the effects of temporary immersion in water).
For specifications such as storage temperature, measurement accuracy, etc., refer to the Instructions for Use that accompanies the equipment used in connection with the accessory.
Wavelength emitted by the sensor: 660 nm and 940 nm in the corresponding Gebrauchsanweisung, die des vom Zähler verwendeten Gerätes beigefügt ist.
The maximum photoic output consumption of the sensor is less than 18 mW.
The information about the wavelength range and maximum photoic output consumption can be especially useful to clinicians ( for example, when photodynamic therapy is performed).

	Authorized representative in the European Community		Date of manufacture
	Not made with natural rubber latex		Batch code
	Ref to instruction manual/booklet		Manufacturer
	Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive).		Protection against the effects of temporary immersion in water

**Rx Only**
**Federal law (USA)** restricts this device to sale by or on the order of a physician.

	Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive).		Protection against the effects of temporary immersion in water
<b>Rx Only</b>	<b>Federal law (USA)</b> restricts this device to sale by or on the order of a physician.		<b>Caution</b>
	Comply with the requirements of Directive 2012/19/EU Waste Electrical & Electronic Equipment		

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